

Anti-POLDIP2 Antibody Picoband® HRP Conjugated

Catalog Number: A09581-4-HRP

About POLDIP2

Polymerase delta-interacting protein 2 also known as Polymerase delta-interacting protein of 38 kDa (PDIP38) is encoded by the POLDIP2 gene in humans. This gene encodes a protein that interacts with the DNA polymerase delta p50 subunit, as well as with proliferating cell nuclear antigen. The encoded protein maybe play a role in the ability of the replication fork to bypass DNA lesions. Alternative splicing results in multiple transcript variants.

Overview

| | |
|----------------------|---|
| Product Name | Anti-POLDIP2 Antibody Picoband® HRP Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ . |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. |
| Host | Rabbit |
| Uniprot ID | Q9Y2S7 |

Technical Details

| | |
|---------------------|---|
| Immunogen | E.coli-derived human POLDIP2 recombinant protein (Position: Q81-Q324). Human POLDIP2 shares 99.2% amino acid (aa) sequence identity with mouse POLDIP2. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | HRP |
| Suggested Dilutions | The intended application should be selected according to the customer's experimental requirements. |

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-POLDIP2 Antibody - HRP

For Research Use Only. Not for use in diagnostic procedures.