

Anti-PEF1 Antibody Picoband® Biotin Conjugated

Catalog Number: A09591-2-Biotin

About PEF1

Peflin is a protein that in humans is encoded by the PEF1 gene. This gene encodes a calcium-binding protein belonging to the penta-EF-hand protein family. The encoded protein has been shown to form a heterodimer with the programmed cell death 6 gene product and may modulate its function in Ca (2+) signaling. Alternative splicing results in multiple transcript variants and a pseudogene has been identified on chromosome 1.

Overview

Product Name	Anti-PEF1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9UBV8

Technical Details

Immunogen	E.coli-derived human PEF1 recombinant protein (Position: R170-L284).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PEF1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.