

## Anti-DGK-Delta (S66) DGKD Antibody

Catalog Number: A09644S66


### Overview

Product Name	Anti-DGK-Delta (S66) DGKD Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DGK-Delta (S66) DGKD Antibody catalog # A09644S66. Tested in WB,IHC,IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q16760


### Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 34-88 of Human DGK-delta.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 IF: 1:50-1:200

## Anti-DGK-Delta (S66) DGKD Antibody (A09644S66) Images



Western blot (WB) analysis of DGK-delta (S66) pAb at 1:500 dilution Lane1:A549 whole cell lysate(40ug) Lane2:Hela whole cell lysate(40ug) Lane3:A2780 whole cell lysate(40ug) Lane4:PC12 whole cell lysate(40ug) Lane5:3T3-L1 whole cell lysate(40ug)



Immunohistochemistry (IHC) analyzes of DGK-delta (S66) pAb in paraffin-embedded human Rectum carcinoma tissue at 1:50, showing cytoplasmic staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DGK-Delta (S66) DGKD Antibody

For Research Use Only. Not for use in diagnostic procedures.