

Anti-FEZF1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A09653-1-Fluoro488

About FEZF1

FEZ family zinc finger 1 is a protein that in humans is encoded by the FEZF1 gene. This gene encodes a transcriptional repressor that belongs to the zinc finger double domain protein family. The encoded protein is thought to play a role in the embryonic migration of gonadotropin-releasing hormone neurons into the brain. Mutations in this gene are associated with hypogonadotropic hypogonadism-22 with anosmia. Alternative splicing results in multiple transcript variants.

Overview

Product Name	Anti-FEZF1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	A0PJY2

Technical Details

Immunogen	E.coli-derived human FEZF1 recombinant protein (Position: A45-K258).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FEZF1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.