

Anti-DCLK2 Antibody

Catalog Number: A09671-2

About DCLK2

DCLK2 belongs to the protein kinase superfamily, contains two N-terminal doublecortin domains, a C-terminal serine/threonine protein kinase domain, and a serine/proline-rich domain in between the doublecortin and the protein kinase domains. DCLK2 is expressed in the brain, heart and eyes (PMID: 18075264).

Overview

Product Name	Anti-DCLK2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-DCLK2 Antibody catalog # A09671-2. Tested in WB, ICC/IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8N568

Technical Details

Immunogen	E.coli-derived human DCLK2 recombinant protein (Position: Q8-D68).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-DCLK2 Antibody

For Research Use Only. Not for use in diagnostic procedures.