

## Anti-BEX3 Antibody Picoband® Cy3 Conjugated

Catalog Number: A09675-2-Cy3

### About BEX3

Protein BEX3 is a protein that in humans is encoded by the NGFRAP1 gene. NGFRAP1 (nerve growth factor receptor associated protein 1), also called BEX3 (brain-expressed X-linked protein 3), NADE (p75NTR-associated death executor) or FAP-1 (Fas-associated phosphatase-1), is a 24 kDa member of the BEX protein family. A nuclear export signal (aa 77-87) allows NGFRAP1 to export bound p75NTR. It is widely expressed, but highest in ovarian granulosa cells, testis, prostate, seminal vesicle and liver. Except for a 13-19 amino acid (aa) rodent Arg/His region at aa 35 that is not found in humans, human NGFRAP1 shares 94% aa identity with mouse and rat NGFRAP1. In humans, a 101 aa isoform with an alternative start site at aa 11 has been described.

### Overview

Product Name	Anti-BEX3 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na2HPO4, 0.02% NaN3.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q00994

### Technical Details

Immunogen	E.coli-derived human BEX3 recombinant protein (Position: Q39-H103).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm

Suggested Dilutions

Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BEX3 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.