

Anti-UBXN1 Antibody Picoband® Biotin Conjugated

Catalog Number: A09683-1-Biotin

About UBXN1

Enables several functions, including enzyme binding activity; polyubiquitin modification-dependent protein binding activity; and proteasome regulatory particle binding activity. Involved in negative regulation of protein metabolic process. Located in cytosol; endoplasmic reticulum; and nucleoplasm. Part of VCP-NPL4-UFD1 AAA ATPase complex.

Overview

Product Name	Anti-UBXN1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q04323

Technical Details

Immunogen	E.coli-derived human UBXN1 recombinant protein (Position: M1-S286). Human UBXN1 shares 92.7% and 92.3% amino acid (aa) sequence identity with mouse and rat UBXN1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBXN1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.