

Anti-UBXN1 Antibody Picoband® Cy3 Conjugated

Catalog Number: A09683-1-Cy3

About UBXN1

Enables several functions, including enzyme binding activity; polyubiquitin modification-dependent protein binding activity; and proteasome regulatory particle binding activity. Involved in negative regulation of protein metabolic process. Located in cytosol; endoplasmic reticulum; and nucleoplasm. Part of VCP-NPL4-UFD1 AAA ATPase complex.

Overview

| | |
|----------------------|--|
| Product Name | Anti-UBXN1 Antibody Picoband® Cy3 Conjugated |
| Reactive Species | Human, Mouse, Rat |
| Application | Flow Cytometry |
| Clonality | Polyclonal |
| Formulation | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light. |
| Host | Rabbit |
| Uniprot ID | Q04323 |

Technical Details

| | |
|---------------------|--|
| Immunogen | E.coli-derived human UBXN1 recombinant protein (Position: M1-S286). Human UBXN1 shares 92.7% and 92.3% amino acid (aa) sequence identity with mouse and rat UBXN1, respectively. |
| Form | Liquid |
| Concentration | 0.5 mg/mL |
| Purification | Immunogen affinity purified. |
| Conjugate | Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users. |

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-UBXN1 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.