

## Anti-UBXN1 Antibody Picoband® PE Conjugated

Catalog Number: A09683-1-PE

### About UBXN1

Enables several functions, including enzyme binding activity; polyubiquitin modification-dependent protein binding activity; and proteasome regulatory particle binding activity. Involved in negative regulation of protein metabolic process. Located in cytosol; endoplasmic reticulum; and nucleoplasm. Part of VCP-NPL4-UFD1 AAA ATPase complex.

### Overview

Product Name	Anti-UBXN1 Antibody Picoband® PE Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, IP, IHC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q04323

### Technical Details

Immunogen	E.coli-derived human UBXN1 recombinant protein (Position: M1-S286). Human UBXN1 shares 92.7% and 92.3% amino acid (aa) sequence identity with mouse and rat UBXN1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	PE Excitation Wavelength: 566 nm Emission Wavelength: 574 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-UBXN1 Antibody - PE

For Research Use Only. Not for use in diagnostic procedures.