

Anti-PHF12 Antibody Picoband® APC Conjugated

Catalog Number: A09690-1-APC

About PHF12

Enables phosphatidylinositol binding activity and transcription corepressor activity. Involved in negative regulation of DNA-templated transcription. Acts upstream of or within negative regulation of transcription by RNA polymerase II. Located in nucleoplasm. Part of Sin3-type complex and transcription repressor complex.

Overview

Product Name	Anti-PHF12 Antibody Picoband® APC Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q96QT6

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PHF12. Human PHF12 shares 100% amino acid (aa) sequence identity with mouse PHF12.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	APC Excitation Wavelength: 633-647 nm Emission Wavelength: 660 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PHF12 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.