

Anti-PHF12 Antibody Picoband® Biotin Conjugated

Catalog Number: A09690-1-Biotin

About PHF12

Enables phosphatidylinositol binding activity and transcription corepressor activity. Involved in negative regulation of DNA-templated transcription. Acts upstream of or within negative regulation of transcription by RNA polymerase II. Located in nucleoplasm. Part of Sin3-type complex and transcription repressor complex.

Overview

Product Name	Anti-PHF12 Antibody Picoband® Biotin Conjugated
Reactive Species	Human, Mouse, Rat
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q96QT6

Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PHF12. Human PHF12 shares 100% amino acid (aa) sequence identity with mouse PHF12.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PHF12 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.