

Anti-CACTIN Antibody

Catalog Number: A09699-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-CACTIN Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-CACTIN Antibody catalog # A09699-1. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat. |
| Application | ELISA, WB |
| Clonality | Polyclonal |
| Formulation | PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q8WUQ7 |

Technical Details

| | |
|---------------------|---|
| Immunogen | The antiserum was produced against synthesized peptide derived from the Internal region of human CACTIN. AA range:141-190 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | Immunogen affinity purified |
| Suggested Dilutions | WB 1:500-2000 ELISA 1:10000-20000 |

Anti-CACTIN Antibody (A09699-1) Images



Western Blot analysis of hepg2 cells using Antibody diluted at 1000. Secondary antibody was diluted at 1:20000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-CACTIN Antibody

For Research Use Only. Not for use in diagnostic procedures.