

Anti-Claudin-8 CLDN8 Antibody

Catalog Number: A09846


Overview

Product Name	Anti-Claudin-8 CLDN8 Antibody
Reactive Species	Mouse, Rat
Description	Boster Bio Anti-Claudin-8 CLDN8 Antibody catalog # A09846. Tested in WB,IHC,ICC/IF applications. This antibody reacts with Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P56750


Technical Details

Immunogen	A synthetic peptide of human CLDN17(NP_036263.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:50-1:100 ICC/IF: 1:50-1:100

Anti-Claudin-8 CLDN8 Antibody (A09846) Images



Western blot analysis of extracts of various cell lines, using CLDN17 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 30s.



Immunohistochemistry of paraffin-embedded mouse kidney using CLDN17 antibody at dilution of 1:100. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Claudin-8 CLDN8 Antibody

For Research Use Only. Not for use in diagnostic procedures.