

Anti-Histone H1 (Mono Methyl Lys25) HIST1H1A Antibody

Catalog Number: A09849-1

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Overview

Product Name	Anti-Histone H1 (Mono Methyl Lys25) HIST1H1A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Histone H1 (Mono Methyl Lys25) HIST1H1A Antibody catalog # A09849-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q02539


Technical Details

Immunogen	Synthetic Peptide of Histone H1 (Mono Methyl Lys25)
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit ascites by affinity-chromatography using specific immunogen.

Suggested Dilutions

WB 1:500-1000

Anti-Histone H1 (Mono Methyl Lys25) HIST1H1A Antibody (A09849-1) Images



Western blotting validation for Anti-Histone H1 (Mono Methyl Lys25) HIST1H1A Antibody A09849-1 Western Blot (WB) analysis of 1) HeLa

[Anti-Histone H1 \(Mono Methyl Lys25\) HIST1H1A Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.