

Anti-RBAK (W22) Antibody

Catalog Number: A09876-1


Overview

| | |
|----------------------|--|
| Product Name | Anti-RBAK (W22) Antibody |
| Reactive Species | Human, Mouse, Rat |
| Description | Boster Bio Anti-RBAK (W22) Antibody catalog # A09876-1. Tested in WB applications. This antibody reacts with Human,Mouse,Rat. |
| Application | WB |
| Clonality | Polyclonal |
| Formulation | Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2 |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9NYW8 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Synthetic peptide, corresponding to amino acids 10-70 of Human RBAK. |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE). |
| Suggested Dilutions | WB: 1:500-1:1000 |

Anti-RBAK (W22) Antibody (A09876-1) Images



Western blot (WB) analysis of RBAK (W22) polyclonal antibody at 1:500 dilution Lane1:L02 whole cell lysate(40ug) Lane2:H1792 whole cell lysate(40ug) Lane3:HEK293T whole cell lysate(40ug) Lane4:PC12 whole cell lysate(40ug) Lane5:AML-12 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RBAK (W22) Antibody

For Research Use Only. Not for use in diagnostic procedures.