

Anti-ApoO Monoclonal Antibody

Catalog Number: A09879


Overview

Product Name	Anti-ApoO Monoclonal Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-ApoO Monoclonal Antibody catalog # A09879. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, IF, IHC, WB
Clonality	Monoclonal 2F1
Formulation	Ascitic fluid containing 0.03% sodium azide.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Mouse
Uniprot ID	Q9BUR5


Technical Details

Immunogen	Purified recombinant fragment of ApoO expressed in E. Coli.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Ascitic fluid
Suggested Dilutions	WB 1:500-1:2000 IHC 1:200-1:1000 ELISA 1:10000 IF 1:50-200

Anti-ApoO Monoclonal Antibody (A09879) Images



Western Blot analysis using ApoO Monoclonal Antibody against HepG2 (1) and 3T3L1 (2) cell lysate



Immunohistochemistry analysis of paraffin-embedded human Testis tissues with AEC staining using ApoO Monoclonal Antibody.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ApoO Monoclonal Antibody

For Research Use Only. Not for use in diagnostic procedures.