

## Anti-PI-16/PI16 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A09887-1-Fluoro550

### About PI16

Peptidase inhibitor 16 is a protein that in humans is encoded by the PI16 gene. PI16 (Peptidase Inhibitor 16) is a variable molecular weight (MW) member of the CRISP family of proteins. It is expressed by cardiomyocytes and serves as an autocrine negative growth regulator.

### Overview

Product Name	Anti-PI-16/PI16 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Monkey, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q6UXB8

### Technical Details

Immunogen	E.coli-derived human PI-16/PI16 recombinant protein (Position: L28-L363).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-PI-16/PI16 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.