

## Anti-LRIG2 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A09895-Fluoro594

### About LRIG2

This gene encodes a transmembrane protein containing leucine-rich repeats and immunoglobulin-like domains. The encoded protein promotes epidermal growth factor signalling, resulting in increased proliferation. Its expression in the cytoplasm of glioma cells is correlated with poor survival. Mutations in this gene can cause urofacial syndrome. Alternative splicing results in multiple transcript variants.

### Overview

Product Name	Anti-LRIG2 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O94898

### Technical Details

Immunogen	E.coli-derived human LRIG2 recombinant protein (Position: E108-R977).
Cross Reactivity	No cross reactivity with other proteins.
Isotype	IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LRIG2 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.