

Anti-TGM4 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A09912-1-Fluoro488

About TGM4

Predicted to enable protein-glutamine gamma-glutamyltransferase activity. Predicted to act upstream of or within mating plug formation. Located in Golgi apparatus.

Overview

Product Name	Anti-TGM4 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P49221

Technical Details

Immunogen	E.coli-derived human TGM4 recombinant protein (Position: M1-K684). Human TGM4 shares 53.3% and 54.2% amino acid (aa) sequence identity with mouse and rat TGM4, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

Anti-TGM4 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.