

## Anti-PYST2/DUSP7 Antibody Picoband® APC Conjugated

Catalog Number: A09919-1-APC

### About DUSP7

Dual-specificity phosphatases (DUSPs) constitute a large heterogeneous subgroup of the type I cysteine-based protein-tyrosine phosphatase superfamily. DUSPs are characterized by their ability to dephosphorylate both tyrosine and serine/threonine residues. DUSP7 belongs to a class of DUSPs, designated MKPs, that dephosphorylate MAPK (mitogen-activated protein kinase) proteins ERK (see MIM 601795), JNK (see MIM 601158), and p38 (see MIM 600289) with specificity distinct from that of individual MKP proteins. MKPs contain a highly conserved C-terminal catalytic domain and an N-terminal Cdc25 (see MIM 116947)-like (CH2) domain. MAPK activation cascades mediate various physiologic processes, including cellular proliferation, apoptosis, differentiation, and stress responses (summary by Patterson et al., 2009 [PubMed 19228121]).

### Overview

|                      |  |
|----------------------|--|
| Product Name         | Anti-PYST2/DUSP7 Antibody Picoband® APC Conjugated   |
| Reactive Species     | Human, Mouse, Rat  |
| Application          | Flow Cytometry   |
| Clonality            | Polyclonal   |
| Formulation          | Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N. |
| Storage Instructions | At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.         |
| Host                 | Rabbit   |
| Uniprot ID           | Q16829   |

### Technical Details

|                     |   |
|---------------------|---|
| Immunogen           | E.coli-derived human PYST2/DUSP7 recombinant protein (Position: G107-T419). |
| Form                | Liquid  |
| Concentration       | 0.5 mg/mL   |
| Purification        | Immunogen affinity purified.  |
| Conjugate           | APC<br>Excitation Wavelength: 633-647 nm<br>Emission Wavelength: 660 nm     |
| Suggested Dilutions | Flow Cytometry, Optimal dilutions should be determined by end users.        |

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-PYST2/DUSP7 Antibody - APC

For Research Use Only. Not for use in diagnostic procedures.