

## Anti-RRP8 Antibody Picoband® Fluoro594 Conjugated

Catalog Number: A09950-2-Fluoro594

### About RRP8

Ribosomal RNA-processing protein 8 is a protein that in humans is encoded by the RRP8 gene. Enables methylated histone binding activity. Involved in several processes, including cellular response to glucose starvation; intrinsic apoptotic signaling pathway by p53 class mediator; and regulation of gene expression. Located in several cellular components, including cytosol; nuclear lumen; and rDNA heterochromatin. Part of chromatin silencing complex.

### Overview

Product Name	Anti-RRP8 Antibody Picoband® Fluoro594 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	O43159

### Technical Details

Immunogen	E.coli-derived human RRP8 recombinant protein (Position: K199-R456).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro594 Excitation Wavelength: 593 nm Emission Wavelength: 618 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RRP8 Antibody - Fluoro594

For Research Use Only. Not for use in diagnostic procedures.