

## Anti-FAM96B/CIAO2B Antibody

Catalog Number: A09953

### Overview

Product Name	Anti-FAM96B/CIAO2B Antibody
Reactive Species	Human
Description	Boster Bio Anti-FAM96B/CIAO2B Antibody catalog # A09953. Tested in WB, IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN <sub>3</sub> , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt -20°C as supplied. 6 months 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9Y3D0

### Technical Details

Immunogen	E.coli-derived human FAM96B recombinant protein (Position: M1-S160).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunohistochemistry, 1:50-400

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-FAM96B/CIAO2B Antibody

For Research Use Only. Not for use in diagnostic procedures.