

## Anti-PHF23 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A09961-1-Fluoro488

### About PHF23

Predicted to enable zinc ion binding activity. Involved in negative regulation of autophagosome assembly; negative regulation of autophagosome maturation; and positive regulation of protein ubiquitination. Located in cytosol and nucleoplasm.

### Overview

Product Name	Anti-PHF23 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9BUL5

### Technical Details

Immunogen	A synthetic peptide corresponding to a sequence at the N-terminus of human PHF23. Human PHF23 shares 95.5% amino acid (aa) sequence identity with mouse PHF23.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your



reviews help your fellow scientists make the right decisions. Thank you for your contribution.

**Anti-PHF23 Antibody - Fluoro488**

**For Research Use Only. Not for use in diagnostic procedures.**