

## Anti-EBP2/EBNA1BP2 Antibody Picoband® Cy3 Conjugated

Catalog Number: A09966-1-Cy3

### About EBNA1BP2

Enables RNA binding activity. Predicted to be involved in rRNA processing and ribosomal large subunit biogenesis. Located in chromosome and nucleolus.

### Overview

Product Name	Anti-EBP2/EBNA1BP2 Antibody Picoband® Cy3 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% NaN <sub>3</sub> .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q99848

### Technical Details

Immunogen	E.coli-derived human EBP2/EBNA1BP2 recombinant protein (Position: M1-M233).
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Cy3 Excitation Wavelength: 554 nm Emission Wavelength: 568 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-EBP2/EBNA1BP2 Antibody - Cy3

For Research Use Only. Not for use in diagnostic procedures.