

## Anti-hnRNP DL Antibody

Catalog Number: A09982-2

### Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



## Overview

|                      |                                                                                                                                                                                                                                                                                                                                                                                                           |
|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name         | Anti-hnRNP DL Antibody                                                                                                                                                                                                                                                                                                                                                                                    |
| Reactive Species     | Human, Mouse, Rat                                                                                                                                                                                                                                                                                                                                                                                         |
| Description          | Boster Bio Anti-hnRNP DL Antibody catalog # A09982-2. Tested in ELISA, IHC applications. This antibody reacts with Human, Mouse, Rat.                                                                                                                                                                                                                                                                     |
| Application          | ELISA, IF, IHC                                                                                                                                                                                                                                                                                                                                                                                            |
| Clonality            | Polyclonal                                                                                                                                                                                                                                                                                                                                                                                                |
| Formulation          | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.                                                                                                                                                                                                                                                                |
| Host                 | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                    |
| Uniprot ID           | O14979                                                                                                                                                                                                                                                                                                                                                                                                    |

## Technical Details


|               |                                                                                                                       |
|---------------|-----------------------------------------------------------------------------------------------------------------------|
| Immunogen     | The antiserum was produced against synthesized peptide derived from human HNRPD. AA range:241-290                     |
| Isotype       | IgG                                                                                                                   |
| Form          | Liquid                                                                                                                |
| Concentration | 1 mg/ml                                                                                                               |
| Purification  | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |

Suggested Dilutions

IHC 1:100-1:300  
ELISA 1:20000  
IF 1:50-200

## Anti-hnRNP DL Antibody (A09982-2) Images

---



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using HNRPD<sup>L</sup> Antibody. The picture on the right is blocked with the synthesized peptide.

### [Anti-hnRNP DL Antibody](#)

For Research Use Only. Not for use in diagnostic procedures.