

Anti-SERTAD1 Antibody

Catalog Number: A09983


Overview

Product Name	Anti-SERTAD1 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SERTAD1 Antibody catalog # A09983. Tested in WB applications. This antibody reacts with Human,Mouse,Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UHV2

Technical Details

Immunogen	A synthetic peptide corresponding to residues in Human SERTAD1.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000

Anti-SERTAD1 Antibody (A09983) Images



Western blot (WB) analysis of SERTAD1 polyclonal antibody at 1:500 dilution Lane1:Jurkat cell lysate Lane2:sp2/0 cell lysate Lane3:PC12 cell lysate

1 Publications Citing This Product

1. PubMed ID: 10.1111/j.1440-1681.2009.05181.x, CONTINUOUS FLUOXETINE ADMINISTRATION PREVENTS RECURRENCE OF PULMONARY ARTERIAL HYPERTENSION AND PROLONGS SURVIVAL IN RATS

Visit bosterbio.com/anti-sertad1-antibody-a09983-boster.html to see all 1 publications.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SERTAD1 Antibody

For Research Use Only. Not for use in diagnostic procedures.