

Anti-RRAGD Antibody

Catalog Number: A09990-2

About RRAGD

Ras-related GTP-binding protein D (RRAGD) is a monomeric GTP or GDP binding protein, which plays a critical role in mediating amino acid stimulated mTOR signaling pathway, a key pathway in determination of cell growth and proliferation rate. Recently, RRAGD has been reported involving the progression of ovarian cancer and colon cancer.

Overview

Product Name	Anti-RRAGD Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RRAGD Antibody catalog # A09990-2. Tested in WB, ICC/IF, IP, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IP, IF, ICC, WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NQL2

Technical Details

Immunogen	E.coli-derived human RRAGD recombinant protein (Position: A3-R246).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000 Immunocytochemistry/Immunofluorescence, 1:50-400 ImmunoPrecipitation, 1:250-300 ELISA, 1:100-1000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RRAGD Antibody

For Research Use Only. Not for use in diagnostic procedures.