

Anti-ELAV-like protein 3 ELAVL3 Antibody

Catalog Number: A10061


Overview

Product Name	Anti-ELAV-like protein 3 ELAVL3 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ELAV-like protein 3 ELAVL3 Antibody catalog # A10061. Tested in WB,IHC,ICC/IF applications. This antibody reacts with Human,Mouse,Rat.
Application	IF, IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q14576


Technical Details

Immunogen	A synthetic peptide of human ELAVL3(NP_001411.2).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 ICC/IF: 1:50-1:200


Anti-ELAV-like protein 3 ELAVL3 Antibody (A10061) Images




Western blot analysis of various lysates, using ELAVL3 antibody at 1:550 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 0.8s.



Immunofluorescence analysis of SH-SY5Y cells using ELAVL3 Rabbit pAb at dilution of 1:100. Blue: DAPI for nuclear staining.



Immunohistochemistry of paraffin-embedded rat brain using ELAVL3 antibody at dilution of 1:200. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry of paraffin-embedded mouse brain using ELAVL3 antibody at dilution of 1:200. Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ELAV-like protein 3 ELAVL3 Antibody

For Research Use Only. Not for use in diagnostic procedures.