

Anti-SLC51A Antibody

Catalog Number: A10127-1

About SLC51A

Predicted to enable protein heterodimerization activity; protein homodimerization activity; and transmembrane transporter activity. Involved in bile acid secretion. Located in basolateral plasma membrane. Implicated in progressive familial intrahepatic cholestasis.

Overview

Product Name	Anti-SLC51A Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-SLC51A Antibody catalog # A10127-1. Tested in WB applications. This antibody reacts with Human, Mouse, Rat.
Application	WB
Clonality	Polyclonal
Formulation	500 ug/ml antibody with PBS, 0.02% NaN ₃ , 1 mg stabilizing protein and 50% glycerol This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	12 months from date of receipt at -20°C as supplied. 6 months at 2 to 8°C after reconstitution. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q86UW1

Technical Details

Immunogen	E.coli-derived human SLC51A recombinant protein (Position: P3-A45).
Form	Liquid
Concentration	500 ug/ml
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 1:500-2000

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SLC51A Antibody

For Research Use Only. Not for use in diagnostic procedures.