

Anti-RNF15 (E193) TRIM38 Antibody

Catalog Number: A10133-1


Overview

Product Name	Anti-RNF15 (E193) TRIM38 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-RNF15 (E193) TRIM38 Antibody catalog # A10133-1. Tested in WB,IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	O00635

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 171-220 of Human RNF15.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 IHC: 1:100-1:300

Anti-RNF15 (E193) TRIM38 Antibody (A10133-1) Images



Western blot (WB) analysis of RNF15 (E193) pAb at 1:500 dilution Lane1:U-87MG whole cell lysate(40ug) Lane2:HeLa whole cell lysate(40ug) Lane3:H9C2 whole cell lysate(40ug) Lane4:3T3-L1 whole cell lysate(40ug)

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RNF15 (E193) TRIM38 Antibody

For Research Use Only. Not for use in diagnostic procedures.