

Anti-Coronin-1B CORO1B Antibody

Catalog Number: A10210


Overview

Product Name	Anti-Coronin-1B CORO1B Antibody
Reactive Species	Rat
Description	Boster Bio Anti-Coronin-1B CORO1B Antibody catalog # A10210. Tested in WB applications. This antibody reacts with Rat.
Application	WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BR76

Technical Details

Immunogen	Recombinant fusion protein of human CORO1B(NP_065174.1).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000

Anti-Coronin-1B CORO1B Antibody (A10210) Images



Western blot analysis of extracts of rat kidney, using CORO1B antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 10s.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Coronin-1B CORO1B Antibody

For Research Use Only. Not for use in diagnostic procedures.