

## Anti-GPR172B SLC52A1 Antibody

Catalog Number: A10245


### Overview

Product Name	Anti-GPR172B SLC52A1 Antibody
Reactive Species	Human
Description	Boster Bio Anti-GPR172B SLC52A1 Antibody catalog # A10245. Tested in WB, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9NWF4


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human PEVR2. AA range:235-284
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000


## Anti-GPR172B SLC52A1 Antibody (A10245) Images



Western Blot analysis of 293T cells using GPR172B Polyclonal Antibody diluted at 1:1000



Western blot analysis of lysates from LOVO cells, using PEVR2 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunofluorescence analysis of MCF7 cells, using PEVR2 Antibody. The picture on the right is blocked with the synthesized peptide.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-GPR172B SLC52A1 Antibody

For Research Use Only. Not for use in diagnostic procedures.