

Anti-LMOD3 Antibody Picoband®

Catalog Number: A10356-2

About LMOD3

Leiomodin-3 is a protein that in humans is encoded by the LMOD3 gene. The protein encoded by this gene is a member of the leiomodin family of proteins. This protein contains three actin-binding domains, a tropomyosin domain, a leucine-rich repeat domain, and a Wiskott-Aldrich syndrome protein homology 2 domain (WH2). Localization of this protein to the pointed ends of thin filaments has been observed, and there is evidence that this protein acts as a catalyst of actin nucleation, and is important to the organization of sarcomeric thin filaments in skeletal muscles. Mutations in this gene have been associated as one cause of Nemaline myopathy, as other genes have also been linked to this disorder. Nemaline myopathy is a disorder characterized by nonprogressive generalized muscle weakness and protein inclusions (nemaline bodies) in skeletal myofibers. Patients with mutations in this gene often present with a severe congenital form of the disorder.

Overview


Product Name	Anti-LMOD3 Antibody Picoband®
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-LMOD3 Antibody Picoband® catalog # A10356-2. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na ₂ HPO ₄ .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q0VAK6

Technical Details

Immunogen	E.coli-derived human LMOD3 recombinant protein (Position: D45-E558). Human LMOD3 shares 72.2% amino acid (aa) sequence identity with mouse LMOD3.
Recommended Detection Systems	Boster recommends Enhanced Chemiluminescent Kit with anti-Rabbit IgG (EK1002) for Western blot.
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.

Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Mouse, Rat ELISA, 0.1-0.5 ug/ml, -

Anti-LMOD3 Antibody Picoband® (A10356-2) Images



Western blot analysis of LMOD3 using anti-LMOD3 antibody (A10356-2). Electrophoresis was performed on a 5-20% SDS-PAGE gel at 70V (Stacking gel) / 90V (Resolving gel) for 2-3 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: rat skeletal muscle tissue lysates, Lane 2: mouse heart tissue lysates, Lane 3: mouse skeletal muscle tissue lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-LMOD3 antigen affinity purified polyclonal antibody (Catalog # A10356-2) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an Enhanced Chemiluminescent detection (ECL) kit (Catalog # EK1002) with Tanon 5200 system. A specific band was detected for LMOD3 at approximately 70 kDa. The expected band size for LMOD3 is at 65 kDa.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-LMOD3 Antibody

For Research Use Only. Not for use in diagnostic procedures.