

## Anti-17 beta-HSD11 HSD17B11 Antibody

Catalog Number: A10397


### Overview

Product Name	Anti-17 beta-HSD11 HSD17B11 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-17 beta-HSD11 HSD17B11 Antibody catalog # A10397. Tested in WB, IF, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, IF, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q8NBQ5


### Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human DHRS8. AA range:71-120
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:5000 IF 1:100-300


## Anti-17 beta-HSD11 HSD17B11 Antibody (A10397) Images




Western Blot analysis of various cells using 17beta-HSD11 Polyclonal Antibody



Western blot analysis of lysates from 293 cells, using DHRS8 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from RAW264.7 cells using DHRS8 antibody.



Immunofluorescence analysis of A549. 1, primary Antibody was diluted at 1:200 (4°C overnight). 2, Goat Anti Rabbit IgG (H&L) - Alexa Fluor 488 Secondary antibody was diluted at 1:1000 (room temperature, 50min). 3, Picture B: DAPI (blue) 10min.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-17 beta-HSD11 HSD17B11 Antibody

For Research Use Only. Not for use in diagnostic procedures.