

## Anti-PIP4K2C Antibody Picoband®

Catalog Number: A10431-1

### About PIP4K2C

This locus represents naturally occurring read-through transcription between the neighboring PMF1 (polyamine-modulated factor 1) and BGLAP (bone gamma-carboxyglutamate Gla protein) genes on chromosome 1. Alternative splicing results in multiple transcript variants encoding isoforms that share sequence identity with the upstream gene product, but they contain distinct C-termini due to frameshifts versus the downstream gene coding sequence.

### Overview


Product Name	Anti-PIP4K2C Antibody Picoband®
Reactive Species	Human
Description	Boster Bio Anti-PIP4K2C Antibody Picoband® catalog # A10431-1. Tested in WB, ICC, IF applications. This antibody reacts with Human. The brand Picoband indicates this is a premium antibody that guarantees superior quality, high affinity, and strong signals with minimal background in Western blot applications. Only our best-performing antibodies are designated as Picoband, ensuring unmatched performance.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Each vial contains 4 mg Trehalose, 0.9 mg NaCl, 0.2 mg Na <sub>2</sub> HPO <sub>4</sub> .
Storage Instructions	At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q8TBX8

### Technical Details


Immunogen	A synthetic peptide corresponding to a sequence in the middle region of human PIP4K2C.
Form	Lyophilized
Concentration	Adding 0.2 ml of distilled water will yield a concentration of 500 ug/ml.
Purification	Immunogen affinity purified.
Suggested Dilutions	Western blot, 0.25-0.5 ug/ml, Human Immunocytochemistry/Immunofluorescence, 5 ug/ml, Human Flow Cytometry (Fixed), 1-3 ug/1x10 <sup>6</sup> cells, Human



## Anti-PIP4K2C Antibody Picoband® (A10431-1) Images



Western blot analysis of PIP4K2C using anti-PIP4K2C antibody (A10431-1). Electrophoresis was performed on a 10% SDS-PAGE gel at 80V (Stacking gel) / 120V (Resolving gel) for 2 hours. The sample well of each lane was loaded with 30 ug of sample under reducing conditions. Lane 1: human RT4 whole cell lysates, Lane 2: human MCF-7 whole cell lysates, Lane 3: human Caco-2 whole cell lysates, Lane 4: human SIHA whole cell lysates. After electrophoresis, proteins were transferred to a nitrocellulose membrane at 150 mA for 50-90 minutes. Blocked the membrane with 5% non-fat milk/TBS for 1.5 hour at RT. The membrane was incubated with rabbit anti-PIP4K2C antigen affinity purified polyclonal antibody (A10431-1) at 0.5 ug/mL overnight at 4°C, then washed with TBS-0.1%Tween 3 times with 5 minutes each and probed with a goat anti-rabbit IgG-HRP secondary antibody at a dilution of 1:5000 for 1.5 hour at RT. The signal is developed using an ECL Plus Western Blotting Substrate (Catalog # AR1196-200) with Tanon 5200 system. A specific band was detected for PIP4K2C at approximately 47 kDa. The expected band size for PIP4K2C is at 47 kDa.



IF analysis of PIP4K2C using anti-PIP4K2C antibody (A10431-1) and anti-Alpha Tubulin antibody (M03989-3). PIP4K2C was detected in an immunocytochemical section of A549 cells. Enzyme antigen retrieval was performed using IHC enzyme antigen retrieval reagent (AR0022) for 15 mins. The cells were blocked with 10% goat serum. And then incubated with 5 ug/mL rabbit anti-PIP4K2C Antibody (A10431-1) and mouse anti-Alpha Tubulin antibody (M03989-3) overnight at 4°C. Fluoro488 Conjugated Goat Anti-Rabbit IgG (BA1127) and Cy3 Conjugated Goat Anti-Mouse IgG (BA1031) were used as secondary antibody at 1:500 dilution and incubated for 30 minutes at 37°C. Visualize using a fluorescence microscope and filter sets appropriate for the label used.

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



### Anti-PIP4K2C Antibody

For Research Use Only. Not for use in diagnostic procedures.