

## Anti-RNF144b/Ibrdc2 Antibody

Catalog Number: A10490

### Overview

Product Name	Anti-RNF144b/Ibrdc2 Antibody
Reactive Species	Human, Mouse
Description	Boster Bio Anti-RNF144b/Ibrdc2 Antibody catalog # A10490. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q7Z419


### Technical Details

Immunogen	Synthesized peptide derived from the Internal region of human p53RFP.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:40000

## Anti-RNF144b/Ibrdc2 Antibody (A10490) Images



Western Blot analysis of various cells using p53RFP Polyclonal Antibody



Western Blot analysis of HeLa cells using p53RFP Polyclonal Antibody

### Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-RNF144b/Ibrdc2 Antibody

For Research Use Only. Not for use in diagnostic procedures.