

## Anti-Cytochrome P450 2A7 CYP2A7 Antibody

Catalog Number: A10507


### Overview

Product Name	Anti-Cytochrome P450 2A7 CYP2A7 Antibody
Reactive Species	Human
Description	Boster Bio Anti-Cytochrome P450 2A7 CYP2A7 Antibody catalog # A10507. Tested in WB, IHC, IF, ELISA applications. This antibody reacts with Human.
Application	ELISA, IF, IHC, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P20853


### Technical Details

Immunogen	Synthesized peptide derived from the C-terminal region of human CYP2A7.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:40000 IF 1:50-200

## Anti-Cytochrome P450 2A7 CYP2A7 Antibody (A10507) Images



Western Blot analysis of 22RV1 cells using CYP2A7 Polyclonal Antibody diluted at 1:2000



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200 (4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200 (room temperature, 45min).

### Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome P450 2A7 CYP2A7 Antibody

For Research Use Only. Not for use in diagnostic procedures.