

Anti-Cytochrome P450 2A7 CYP2A7 Antibody

Catalog Number: A10507-1


Overview

Product Name	Anti-Cytochrome P450 2A7 CYP2A7 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Cytochrome P450 2A7 CYP2A7 Antibody catalog # A10507-1. Tested in WB,IHC,ICC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P20853

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 435-480 of Human CYP2A7.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200 ICC: 1:50-1:200

Anti-Cytochrome P450 2A7 CYP2A7 Antibody (A10507-1) Images



Western blot (WB) analysis of CYP2A7 polyclonal antibody at 1:500 dilution Lane1:Hela cell lysate Lane2:sp2/0 cell lysate Lane3:PC12 cell lysate

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Cytochrome P450 2A7 CYP2A7 Antibody

For Research Use Only. Not for use in diagnostic procedures.