

Anti-Inhibin Beta-C (V129) INHBC Antibody

Catalog Number: A10556-2


Overview

Product Name	Anti-Inhibin Beta-C (V129) INHBC Antibody
Reactive Species	Human
Description	Boster Bio Anti-Inhibin Beta-C (V129) INHBC Antibody catalog # A10556-2. Tested in WB,IHC applications. This antibody reacts with Human.
Application	IHC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P55103

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 100-150 of Human Inhibin beta-C.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:1000 IHC: 1:50-1:200

Anti-Inhibin Beta-C (V129) INHBC Antibody (A10556-2) Images



Immunohistochemistry (IHC) analyzes of Inhibin beta-C (V129) pAb in paraffin-embedded human liver carcinoma tissue at 1:50, showing strong cytoplasmic staining. Negative control (the right) Using PBS instead of primary antibody, secondary antibody is Goat Anti-Rabbit IgG-biotin followed by avidin-peroxidase.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-Inhibin Beta-C (V129) INHBC Antibody

For Research Use Only. Not for use in diagnostic procedures.