

Anti-TWIK-2 KCNK6 Antibody

Catalog Number: A10559


Overview

| | |
|----------------------|---|
| Product Name | Anti-TWIK-2 KCNK6 Antibody |
| Reactive Species | Human |
| Description | Boster Bio Anti-TWIK-2 KCNK6 Antibody catalog # A10559. Tested in ELISA, WB applications. This antibody reacts with Human. |
| Application | ELISA, WB |
| Clonality | Polyclonal |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required. |
| Storage Instructions | Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles. |
| Host | Rabbit |
| Uniprot ID | Q9Y257 |

Technical Details

| | |
|---------------------|---|
| Immunogen | Synthesized peptide derived from TWIK-2 . at AA range: 240-320 |
| Isotype | IgG |
| Form | Liquid |
| Concentration | 1 mg/ml |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Suggested Dilutions | WB 1:500-1:2000 ELISA 1:5000 |

Anti-TWIK-2 KCNK6 Antibody (A10559) Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4°over night

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TWIK-2 KCNK6 Antibody

For Research Use Only. Not for use in diagnostic procedures.