

## Anti-TSPYL5 Antibody Picoband® Biotin Conjugated

Catalog Number: A10648-2-Biotin

### About TSPYL5

TSPYL5 is a protein that in humans is encoded by the TSPYL5 gene. Predicted to enable chromatin binding activity and histone binding activity. Involved in several processes, including cellular response to gamma radiation; positive regulation of protein kinase B signaling; and positive regulation of protein ubiquitination. Predicted to be active in chromatin and nucleus.

### Overview

Product Name	Anti-TSPYL5 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Application	WB, IHC, ELISA
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na <sub>2</sub> HPO <sub>4</sub> , 0.02% Na <sub>3</sub> N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q86VY4

### Technical Details

Immunogen	E.coli-derived human TSPYL5 recombinant protein (Position: R24-Q337).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	Western blot, Optimal dilutions should be determined by end users. Immunohistochemistry (Paraffin-embedded Section), Optimal dilutions should be determined by end users. ELISA, Optimal dilutions should be determined by end users.

## Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TSPYL5 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.