

Anti-TSPYL5 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A10648-2-Fluoro488

About TSPYL5

TSPYL5 is a protein that in humans is encoded by the TSPYL5 gene. Predicted to enable chromatin binding activity and histone binding activity. Involved in several processes, including cellular response to gamma radiation; positive regulation of protein kinase B signaling; and positive regulation of protein ubiquitination. Predicted to be active in chromatin and nucleus.

Overview

Product Name	Anti-TSPYL5 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q86VY4

Technical Details

Immunogen	E.coli-derived human TSPYL5 recombinant protein (Position: R24-Q337).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TSPYL5 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.