

Anti-SNRPA1 Antibody Picoband® Fluoro550 Conjugated

Catalog Number: A10655-1-Fluoro550

About SNRPA1

U2 small nuclear ribonucleoprotein A' is a protein that in humans is encoded by the SNRPA1 gene. Enables RNA binding activity. Involved in mRNA splicing, via spliceosome and spermatogenesis. Located in nuclear speck. Part of U2-type catalytic step 2 spliceosome and U2-type pre-catalytic spliceosome. Implicated in connective tissue disease.

Overview

Product Name	Anti-SNRPA1 Antibody Picoband® Fluoro550 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Flow Cytometry
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	P09661

Technical Details

Immunogen	E.coli-derived human SNRPA1 recombinant protein (Position: M1-E233).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro550 Excitation Wavelength: 562 nm Emission Wavelength: 576 nm
Suggested Dilutions	Flow Cytometry, Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SNRPA1 Antibody - Fluoro550

For Research Use Only. Not for use in diagnostic procedures.