

Anti-Synaptotagmin-11 SYT11 Antibody

Catalog Number: A10695-2

Submit a product review to Biocompare.com

[Submit a review of this product](#) to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.




Overview

Product Name	Anti-Synaptotagmin-11 SYT11 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-Synaptotagmin-11 SYT11 Antibody catalog # A10695-2. Tested in WB, ICC/IF applications. This antibody reacts with Human, Mouse, Rat.
Application	IF, ICC, WB
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9BT88


Technical Details

Immunogen	Recombinant fusion protein of human SYT11(NP_689493.3).
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	WB: 1:500-1:2000 ICC/IF: 1:50-1:100


Anti-Synaptotagmin-11 SYT11 Antibody (A10695-2) Images




Western blot analysis of extracts of various cell lines, using SYT11 antibody at 1:670 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 20s.



Immunofluorescence analysis of mouse brain cells using SYT11 Rabbit pAb at dilution of 1:150. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of rat brain cells using SYT11 Rabbit pAb at dilution of 1:150. Blue: DAPI for nuclear staining.



Immunofluorescence analysis of SH-5Y5Y cells using SYT11 Rabbit pAb at dilution of 1:150. Blue: DAPI for nuclear staining.

Anti-Synaptotagmin-11 SYT11 Antibody

For Research Use Only. Not for use in diagnostic procedures.