

Anti-BMP-8 (E90) BMP8B Antibody

Catalog Number: A10747-1


Overview

Product Name	Anti-BMP-8 (E90) BMP8B Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-BMP-8 (E90) BMP8B Antibody catalog # A10747-1. Tested in IHC applications. This antibody reacts with Human,Mouse,Rat.
Application	IHC
Clonality	Polyclonal
Formulation	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH7.2
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	P34820/Q7Z5Y6

Technical Details

Immunogen	Synthetic peptide, corresponding to amino acids 70-120 of Human BMP-8.
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
Suggested Dilutions	IHC: 1:50-1:200

Anti-BMP-8 (E90) BMP8B Antibody (A10747-1) Images



Immunohistochemistry (IHC) analyzes of BMP-8 (E90) pAb in paraffin-embedded human testis tissue.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-BMP-8 (E90) BMP8B Antibody

For Research Use Only. Not for use in diagnostic procedures.