

Anti-MAN2C1 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A10839-Fluoro488

About MAN2C1

MAN2C1 belongs to the glycosyl hydrolase 38 family.

Overview

Product Name	Anti-MAN2C1 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human, Mouse, Rat
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, IF, ICC, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9NTJ4

Technical Details

Immunogen	E.coli-derived human MAN2C1 recombinant protein (Position: A176-Q1027). Human MAN2C1 shares 90.5% and 87.3% amino acid (aa) sequence identity with mouse and rat MAN2C1, respectively.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-MAN2C1 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.