

Anti-SHROOM4 Antibody Picoband® Biotin Conjugated

Catalog Number: A10866-Biotin

About SHROOM4

This gene encodes a member of the APX/Shroom family, which contain an N-terminal PDZ domain and a C-terminal ASD2 motif. The encoded protein may play a role in cytoskeletal architecture. Mutations in this gene have been linked to the X-linked Stocco dos Santos syndrome characterized by cognitive disabilities. Alternatively spliced transcript variants have been described.

Overview

Product Name	Anti-SHROOM4 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% Na ₃ N.
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9ULL8

Technical Details

Immunogen	E.coli-derived human SHROOM4 recombinant protein (Position: Y262-H928).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SHROOM4 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.