

Anti-SHROOM4 Antibody Picoband® Fluoro488 Conjugated

Catalog Number: A10866-Fluoro488

About SHROOM4

This gene encodes a member of the APX/Shroom family, which contain an N-terminal PDZ domain and a C-terminal ASD2 motif. The encoded protein may play a role in cytoskeletal architecture. Mutations in this gene have been linked to the X-linked Stocco dos Santos syndrome characterized by cognitive disabilities. Alternatively spliced transcript variants have been described.

Overview

Product Name	Anti-SHROOM4 Antibody Picoband® Fluoro488 Conjugated
Reactive Species	Human
Application	Recommended applications are based on the parent unconjugated antibody (ELISA, Flow Cytometry, WB). Customers may select suitable applications according to their experimental needs.
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Host	Rabbit
Uniprot ID	Q9ULL8

Technical Details

Immunogen	E.coli-derived human SHROOM4 recombinant protein (Position: Y262-H928).
Cross Reactivity	No cross-reactivity with other proteins.
Isotype	Rabbit IgG
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Fluoro488 Excitation Wavelength: 488 nm Emission Wavelength: 515-545 nm
Suggested Dilutions	Optimal dilutions should be determined by end users.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-SHROOM4 Antibody - Fluoro488

For Research Use Only. Not for use in diagnostic procedures.