

Anti-TUFT1 Antibody Picoband® Biotin Conjugated

Catalog Number: A10868-1-Biotin

About TUFT1

Tuftelin is an acidic protein that is thought to play a role in dental enamel mineralization and is implicated in caries susceptibility. It is also thought to be involved with adaptation to hypoxia, mesenchymal stem cell function, and neurotrophin nerve growth factor mediated neuronal differentiation.

Overview

Product Name	Anti-TUFT1 Antibody Picoband® Biotin Conjugated
Reactive Species	Human
Clonality	Polyclonal
Formulation	Each vial contains 50% glycerol, 0.9% NaCl, 0.2% Na ₂ HPO ₄ , 0.02% NaN ₃ .
Storage Instructions	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing.
Host	Rabbit
Uniprot ID	Q9NNX1

Technical Details

Immunogen	E.coli-derived human TUFT1 recombinant protein (Position: T9-E339). Human TUFT1 shares 92.4% amino acid (aa) sequence identity with mouse TUFT1.
Form	Liquid
Concentration	0.5 mg/mL
Purification	Immunogen affinity purified.
Conjugate	Biotin
Suggested Dilutions	The intended application should be selected according to the customer's experimental requirements.

Submit a product review to Biocompare.com

Submit a review of this product to Biocompare.com to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-TUFT1 Antibody - Biotin

For Research Use Only. Not for use in diagnostic procedures.