

Anti-ABCF2 Antibody

Catalog Number: A10873


Overview

Product Name	Anti-ABCF2 Antibody
Reactive Species	Human, Mouse, Rat
Description	Boster Bio Anti-ABCF2 Antibody catalog # A10873. Tested in WB, ELISA applications. This antibody reacts with Human, Mouse, Rat.
Application	ELISA, WB
Clonality	Polyclonal
Formulation	Liquid in PBS containing 50% glycerol, 0.5% stabilizing protein and 0.02% sodium azide. This antibody is supplied in a stabilized formulation. Compatibility with conjugation reactions depends on the chemistry of the conjugation method used. For conjugation methods that are not compatible with the stabilizing components present in this formulation, a carrier-free antibody format is required.
Storage Instructions	Store at -20°C for one year. For short term storage and frequent use, store at 4°C for up to one month. Avoid repeated freeze-thaw cycles.
Host	Rabbit
Uniprot ID	Q9UG63


Technical Details

Immunogen	The antiserum was produced against synthesized peptide derived from human ABCF2. AA range:171-220
Isotype	IgG
Form	Liquid
Concentration	1 mg/ml
Purification	Immunogen affinity purified
Suggested Dilutions	WB 1:500-1:2000 ELISA 1:5000


Anti-ABCF2 Antibody (A10873) Images




Western Blot analysis of various cells using ABCF2 Polyclonal Antibody



Western Blot analysis of HeLa cells using ABCF2 Polyclonal Antibody



Western blot analysis of ABCF2 Antibody. The lane on the right is blocked with the ABCF2 peptide.



Western blot analysis of the lysates from HepG2 cells using ABCF2 antibody.

Submit a product review to [Biocompare.com](https://www.biocompare.com)

Submit a review of this product to [Biocompare.com](https://www.biocompare.com) to receive a \$20 Amazon.com giftcard! Your reviews help your fellow scientists make the right decisions. Thank you for your contribution.



Anti-ABCF2 Antibody

For Research Use Only. Not for use in diagnostic procedures.